





medicalbees STRAIGHT SAW

Jetzt verfügbar im medical bees Produktprogramm

Die neue Stabsäge basierend auf dem bewährten **beesystem** steht jetzt zur Verfügung.

- ✓ Gerades Design (statt abgewinkelten Design) für paralleles Arbeiten mit der oszillierenden Säge
- ✓ Gleiches Powerpack/ Ladegerät/ Zubehör wie für das beesystem können verwendet werden.
- ✓ Gewohntes Handling
- ✓ stufenlos regulierbare Drehzahl in 2 Bereichen
- ✓ Das bewährte Konzept des beesystem bleibt weiterhin bestehen



Art.Nr. /Art.No.	Grad / degree	Für Sägeblätter/ for sawblades	
F-40-100-00	4°	Stryker6 design	
F-40-101-00	2°	Stryker6 design	
F-40-110-00	4°	Synthes/AO design	
F-40-111-00	2°	Synthes/AO design	

The new straight saw based on the proven **beesystem** is now available

- ✓ Straight design (instead of pistol grip) for parallel working with the oscillating saw
- ✓ The same Powerpack/ charger/ accessory as for the beesystem can be used.
- ✓ Handling remains the same.
- ✓ Continuously adjustable rotation speed in 2 steps
- ✓ The proven concept of the beesystem remains in place